

mylife Loop







mylife Loop powered by CamAPS FX


The hybrid closed loop (HCL) system that runs from your smartphone



1. The product images are for illustrative purposes only. mylife and YpsoPump are registred trademarks of Ypsomed AG. CamAPS is a registered trademark of CamDiab Ltd. Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott.
















Enjoy exercise with 'Ease-off'



- Reduces insulin delivery and raises PGT for periods of low glucose
- Adaptive, easy to turn off, and will keep you in range
- Immediately activate or pre-programme
 - Run up to 24 hours at a time




Manage the unexpected with 'Boost'

•



- Increases insulin delivery
- Immediately activate or pre-programme
- Run up to 13 hours at a time
- Adaptive, easy to turn off, and will keep you in range






Make it your own with the Personal Glucose Target



- Programme for better TIR overnight – deal with dawn phenomenon / overnight hypos
- Tight TIR during pregnancy
- Can reduce HbA1c slowly



Take the hassle out of carbs with Add Meal

•



- 'Add Meal' menu makes parties, trips, sleepovers and picky eating easier
- Use 'Slowly Absorbed Meals' for when carb counting might be tricky
 - Easily log hypos treatment too





Testimonials



It is amazing how much treatment has changed since I was diagnosed, I really feel like I've come from the desert into the digital world. To anyone considering the CamAPS FX app **my advice is quite simple, go for it. It really is transformational** and if I can do it anyone can

Roger



"Every day is easier with but I think the best way to evidence it's power is how impacted my yearly trip to Inida. For the first time in my life, my diabetes didn't get in the way of my trip. I was able to forget about it – enjoy the food (fabulous Indian waffles for breakfast), enjoy the heat, ignore the jet lag and be present!"

Phoebe



mylife Loop has changed my life ... from smoothing over blousing errors to sleeping soundly ... being able to trust technology to do things for me is amazing. It's like a weight lifted mentally – there's so much less decision making...it's incredible how much easier it's made coping with T1.

_aura



Clinical research

Developed at the University of Cambridge, the CamAPS FX algorithm is considered by many as the gold standard.

Randomised clinical trials over 15 years have evidenced significant:



For more information on clinical research, please visit the FAQ section on the CamDiab website: <u>https://camdiab.com/faq</u>





SOLUTION

The pump itself is ..



- Small, lightweight and easy to use¹
- lcon based with an intuitive¹ touchscreen
- Waterproof up to 1 meter for 60 minutes
- Powered by a AAA battery ... no wires, no fuss!

The infusion sets ..



- Rotate 360° for flexibility
- Offer a choice of steel or soft cannulas, as well as a variety of tubing lengths
- Can be changed independently of the cartridge

The cartridge options are ..

Ш.	
i i i	
	\square
F	
-"	
T	

- A super easy to change¹ pre-filled cartridge for fast cartridge changes
- Or a self-filled reservoir for insulin of choice

1 Waldenmaier D. et al.: First User Experiences With a Novel Touchscreen-Based Insulin Pump System in Daily Life of Patients With Type 1 Diabetes Experienced in Insulin Pump Therapy. Journal of Diabetes Science and Technology. 2019;13(1):96-102. doi: 10.1177/1932296818785386.

Features & benefits of the mylife YpsoPump



Our service

24/7/365, UK based customer support line with an experienced, confidence and supportive team.



Quick and easy order process over the phone or online.

(alla)	

Access to our clinical education team to ensure best outcomes from mylife Loop.



Exceptional customer service from day one with a small, friendly team.





Download to demo

mylife CamAPS FX demo



mylife YpsoPump Explorer app 3D simulator











More **Freedom.** More **Confidence.** With **mylife.**

